



ProPulse

VersaTron®



GENERAL RECOMMENDATIONS

500-1000 pulses per treatment
1-3 treatments;
delivered 2-3 weeks apart

Other:

Non-union/delayed healing fracture:

5mm or 20mm
850-1000 pulses
E6

TPLO:

5mm or 20mm
850-1000 pulses
1-2 treatments, 1st treatment to be administered at time of surgery

Superficial Wounds

5mm
8 pulses/cm²; min 500
E2

These guidelines are ONLY applicable to PulseVet shock wave; other shock wave devices have a substantially different output.

All guidelines are provided for medium-large dogs with average severity / chronic injuries. For more specific recommendations, please refer to the 'ProPulse Small Animal Application Guide'. The following modifications can be made at the practitioner's discretion:

- For smaller dogs where less surface area needs to be covered, the number of pulses may be reduced slightly and for very large dogs, increased slightly.
- For dogs with acute soft tissue injuries, the energy setting may be reduced slightly.
- For very severe OA and tendonitis cases, a higher number of pulses and energy setting should be used,